




LATERAL
VETERINARY

Warming Mattress + Return Electrode

U501V

19.5 in x 32 in (49.5 cm x 82 cm)³
2.0 lb (0.9 kg)



Easy to Use

- Compatible with most electrosurgical generators⁴
- Use with one or two electrosurgical generators
- Lightweight for quick & easy set up and cleaning
- One product, two functions
- Requires no skin preparation, such as shaving

Improved Safety

- Eliminates risk of burns from loose adhesive pads
- No adhesive - ideal for patients with delicate skin
- Capacitive-coupling return electrodes have been used for 20+ years in over 100 million procedures⁵

Patient Warming

- Air-free, safe, conductive heat
- Better temperature outcomes when combining a Mattress and Blanket

Cost Savings

- Reusable for up to 2 years
- Saves money and time compared to single-use pads or other reusable return electrodes
- Reduces SKUs and ordering hassles

Environmentally Friendly

- Reduces waste from single-use adhesive pads and single-use warming devices
- Lowers pollution from shipping and transport
- Meets green initiatives⁵

1. FDA cleared, 510(k) K230866, 510(k) K232627

2. Megadyne and Mega Soft are trademarks of J&J Ethicon

3. Heated dimensions

4. Adapter may be required

5. jnjmedtech.com/en-US/product/megadyne-mega-soft-patient-return-electrode-pad



Warming Mattress + Return Electrode Cables

A112 (VET)

Warming Mattress Cable


16ft (5 m)



A136 (VET)

Return Electrode Cable


13 ft (4 m)



A137 (VET)

Return Electrode, Dual Generator “Y” Cable

11.5 ft, 4.9 ft (Y Lengths)
(3.5 m, 1.5 m)



A138 (VET)

REM Enabled Return Electrode Adapter

16.5 in (42.0 cm)



LATERAL
VETERINARY

1300 144 734

tristanl@latvet.com | www.latvet.com